

## DATASHEET

---

Ambu® Open Cuff Silicone Face Masks are designed for use with manual and automatic resuscitators and ventilators.

The masks are made of soft silicone and shaped for a comfortable fit for both adults and children.

# Ambu® Open Cuff Silicone Face Mask



## KEY BENEFITS

---

Transparent for visual check of bleeding, vomitus, and spontaneous breathing

Available in 6 different sizes

All models can be autoclaved repeatedly at 134 °C

# Ambu

## MATERIALS

Product	Material
Mask	Silicone rubber

## SIZE

Type	Weight
0 Infant	13 g
1 Infant	26 g
2 Child	41 g
3 Child	54 g
4 Adult/Child	76 g
5 Adult	102 g

## ORDERING SPECIFICATIONS

Item no.	Type
000 251 001	Ambu Open Cuff Silicone Face Mask size 0 Infant
000 251 002	Ambu Open Cuff Silicone Face Mask size 1 Infant
000 251 003	Ambu Open Cuff Silicone Face Mask size 2 Child
000 251 004	Ambu Open Cuff Silicone Face Mask size 3 Child
000 251 005	Ambu Open Cuff Silicone Face Mask size 4 Adult/child
000 251 006	Ambu Open Cuff Silicone Face Mask size 5 Adult



Ambu Open Cuff Silicone Face Mask sizes 0, 1 and 2



Ambu Open Cuff Silicone Face Mask sizes 3, 4 and 5

# Ambu

### Ambu A/S

Baltorpbakken 13  
2750 Ballerup  
Denmark  
T +45 72 25 20 00  
ambu.com

CE 0434 US: Rx Only